


Deluxe Safety Belt Standard

193

- Used for rehabilitation to help maintain patient's mobility
- Often used in combination with other transfer accessories
- Can be used to assist with exercise and mobility
- Non-slip fabric on inner belt to help keep in position and prevent slipping
- Wash temp 71°C / 159°F
- Do not tumble dry



Product Codes & Specifications:



	Polyester and Nylon	A	B	C
Small to Medium	193S2M	87 - 107cm	90cm	15cm
Standard	193SND	97 - 117cm	100cm	15cm
Medium to Large	193M2L	112 - 132cm	106cm	16cm

