


## Back and Leg Raiser

### 201

- Designed to assist carer with supporting a patient
- Helps to reduce handling
- Can be used across shoulders to gain access to insert a sling
- Can be used to assist a patient into a sitting position when in bed or chair
- Can be used either on one leg or both together to raise or reposition and assist with wound or care treatment
- Option of both inner and outer positioning handles
- This is a manual device and current moving and handling guidelines should be followed
- Manufactured in a quilted, polyester fabric with an inner soft lining for additional comfort
- Wash temp 71°C / 159°F
- Tumble dry low heat (max 50°C)



### Product Codes & Specifications:

			
	Polyester	A	B
ONE SIZE	201	58cm	22cm

