

## **Bariatric Flat Slide Sheets**

- Used (single or in pairs) to help slide the patient up or across a surface
- · Reduces the risk of shearing to patient's skin
- Reduces the risk of injury to both users and caregivers
- Used for assisting the insertion and retrieval of a sling
- Manufactured from a low friction material which provides gliding properties and reduced shearing forces
- Essential part of a caregiver's kit
- Wash temp 71°C / 159°F
- Do not tumble dry

## **B-FLAT**



Colours may vary

## **Product Codes & Specifications:**

