

## Swift® Slider

- A friction-reducing device that remains in-situ to assist with care receiver repositioning and turning in bed
- The Swift® In-Bed Repositioning System features two elements, the Swift® Slider and the UltraSlide® Bottom Sheet, which work together to support repositioning and turning manoeuvres in bed
- The breakthrough WhooTex® fabric ensures a soft, comfortable feel against the skin
- The low friction backing allows care practitioners to reposition a patient safely and easily whilst the breathable construction helps to avoid heat build-up
- The Swift® Slider's thin profile will not impair the performance of alternating air pressure reducing mattresses and can therefore remain on the bed when not in use (tucked in), making it readily available when required
- Swift® Slider with Handles run inset down two lengths and can be used to hook up to a hoist spreader bar
- Offers multidirectional movement
- All slider sizes are suitable for both standard and bariatric application

## SIL-PTD-56/BA



Swift® Slider with Handles

### Product Codes & Specifications:

	 SWL	 WhooTex® Fabric	A	B
Standard	400kg	SIL-PTD-56/BA/Y/UK	186cm	123cm
Long	400kg	SIL-PTD-56/BA/LG/UK	205cm	200cm
Long with handles	200kg	SIL-PTD-56/BA/H/UK	243cm	200cm

