

Bradbury Turner

- Easy to apply
- Enables caregivers to turn a patient onto their side with minimum handling
- Turning can be carried out either side of the bed without re-inserting the sling
- Webbing loops on all four corners of the Bradbury Turner allow the device to be attached to hoist from either side of the bed
- Designed to turn and support the patient on their side
- Clunk-click buckle secures to bed frame providing greater safety
- Wash temp 85°C / 185°F
- Tumble dry low heat (max 50°C)

Product Codes & Specifications:

				
100% Polyester		A	B	C
One Size	CC390	up to 158cm	108cm	46cm

CC390



NB: Colours may vary

