

Walk Assist Jacket

- Designed to assist an ambulant patient with an unstable gait
- The jacket has soft handles positioned to enable the carers to offer some postural support and avoid grabbing clothing / limbs
- Carers can assist from behind or at the sides
- Ideal to assist those using walking sticks / frame
- Easy to apply like a jacket with zip fastening and side zips to allow expansion for outer clothing
- Mesh wash bag supplied for use when laundering jacket
- Wash temp 85°C / 185°F
- Tumble dry low heat (max 50°C)



SIL-WALK



Product Codes & Specifications:



	100% Polyester Mesh	A	B	C	D
XXS	SIL-WALK-XXSM	59 - 77cm	38 - 44cm	41cm	53cm
XS	SIL-WALK-XSM	77 - 87cm	44 - 50cm	45cm	58cm
S	SIL-WALK-SM	87 - 107cm	50 - 57cm	49cm	62cm
M	SIL-WALK-MED	97 - 117cm	57 - 64cm	53cm	72cm
L	SIL-WALK-LRG	112 - 132cm	64 - 69cm	59cm	83cm

