

# Slide Mitt (Pair)

- Reduces the risk of shearing to the patient’s skin
- Used to assist with the insertion of a sling
- Can be used singularly or in pairs
- Manufactured from a low friction material which provides gliding properties and reduced shearing forces
- Wash temp 71°C / 159°F
- Do not tumble dry

# SLDMITT



## Product Codes & Specifications:

|                                                                                     |         |      |      |
|-------------------------------------------------------------------------------------|---------|------|------|
|  |         |      |      |
| 100% Polyester<br>silicone coated                                                   |         | A    | B    |
| One Size                                                                            | SLDMITT | 50cm | 19cm |

