

## Repositioning Sheet

## M921522, M921506

- Used for repositioning and lateral transfers
- Used to help with effective pressure management, to reduce and prevent pressure sores
- Constructed with reinforced tapes from the side of the sling
- Can be used as part of the bed make-up when used in combination with a draw sheet
- Generous size making it suitable for all body dimensions
- Wash temp 95°C / 203°F
- Tumble dry low heat (max 50°C)



### Product Codes & Specifications:

	 SWL 454kg	 SWL 454kg		
	<b>100% Polyester</b>	<b>100% Polyester Mesh</b>	<b>A</b>	<b>B</b>
<b>One Size</b>	<b>M921522</b>	<b>M921506</b>	132cm	247cm

