Universal D-anserclip®



Our Universal D-anserclip® fits all the major hoist models within the UK and European markets, allowing greater interchangeability between varying clip hoists.

The safety "Click-on Click-off" action ensures a physical click is heard on the application of the spreader bar. All moving parts are internal to the clip ensuring greater durability during the lifetime of the patient sling.

Stand Up Sling

Code: D202(SIZE)CP

- Good general stand aid
- Fully padded
- · Hook and loop fastening with clunk click buckles for patient safety
- Ideal for maintaining mobility and for building strength
- Weight bearing ability required

Rigid Stand Aid

Code: **DA-AS112**

- Semi-rigid corset-style stand aid
- · Deep back design providing additional support
- Clunk-click buckle and hook and loop fixing waist support for additional safety
- Suitable for maintaining mobility or rehabilitation
- Padded underarm rolls for extra comfort
- · The patient should have some core stability and weight-bearing ability

High Easy Toileting

Code: **DA-ST3A**

- Easy to apply padded sling
- Built in adjustable head support
- Good access for hygiene purposes
- Good upper body tone required
- Built in waist support with clunk-click fastening
- Full risk assessment required for suitability

High Easy

Code: **HYCP**

- Full body sling with head support
- Divided padded leg sections
- Removable head support stiffeners
- Ideal for general transfers
- Small aperture for hygiene purposes

In-Situ High Easy

Code: **HYAD**

- Available in Silva-Superfine® PLUS, Silva-Superfine with COOLMAX® and Double Layer with soft binding
- Designed to provide full body & head support on transfer
- Adjustable integral head support with removable stiffeners
- Easy to apply with low-friction fabric to leg areas for ease of application
- Once the patient has been transferred, this sling can be left in place, subject to an appropriate risk assessment

All-in-One High Easy

Code: **DA-HYA2**

- Designed to provide full body and head support on transfer
- Apply on the bed or in a supine position
- Suitable for amputee scenarios
- Soft binding as standardNo commode aperture
- Universal clip-fixing

Code. DA-ITIAZ































