


Leg Lifting Support Sling

- Manufactured in Silva-Spacer with COOLMAX® and COOLMAX freshFX® EcoMade
- Offers additional support for patient hygiene and wound treatment
- Helps to reduces patient / carer handling
- Can be used during a transfer or to provide support during a procedure
- Wash temp 95°C / 203°F
- Tumble dry low heat (max 50°C)

B-LEG



Product Codes & Specifications:

 			
Silva-Spacer with COOLMAX® and COOLMAX freshFX® EcoMade			
		A	B
One Size	B-LEG	112cm	21cm

