

Cylindrical Slide Sheets (Taped Edges)

SLDT

- Taped edge to assist with repositioning
- Used to help slide the patient up or across a surface
- Reduces the risk of shearing to patient’s skin
- Reduces the risk of injury to both users and caregivers
- Used for assisting the insertion and retrieval of a sling
- It is an essential part of a caregiver’s kit
- Manufactured from a low friction material which provides gliding properties and reduced shearing forces
- Wash temp 71°C / 159°F
- Do not tumble dry



Product Codes & Specifications:

			
	100% Polyester silicone coated	A	B
Standard	SLDSNDT	122cm	68cm
Long	SLDLT	145cm	68cm
Wide	SLDBT	120cm	100cm

