

In-Situ Repositioning Sheet

SKU: M921539


KEY FEATURES

- Used for repositioning and lateral transfers
- Used to help with effective pressure management, to reduce and prevent pressure sores
- Size L - constructed from 3D Spacer fabric
- Size XL - constructed from 3D Spacer fabric with COOLMAX® and COOLMAX freshFX® EcoMade
- Can be used as part of the bed make-up when used in combination with a draw sheet
- Elastic on corners so that the sheet will stay in place on the mattress
- Generous size making it suitable for all body dimensions



PRODUCT CODES AND SPECIFICATIONS

SIZE	PRODUCT CODE: 3D Spacer	PRODUCT CODE: Spacer with COOLMAX® and COOLMAX FreshFX® EcoMade	SWL	A	B
1	M921539-L		454kg	142cm	198cm
2		M921539-XL	544kg	168cm	198cm



CLEANING INSTRUCTIONS

- Wash temp 85°C / 185°F
- Tumble dry low heat (max 50°C)